

PROVOX Life

LaryTube and HMEs

CE 0413 LaryTube

CE HME

UK
CA HME

LaryTube

Rx
ONLY

LaryTube


HME

42 °C
108 °F
RT
2 °C
36 °F

Provox Life LaryTube Figures


1.


8	12.0mm	27/36/55mm
9	13.5mm	27/36/55mm
10	15.0mm	27/36/55mm
12	17.0mm	27/36/55mm


2.


3.


4a


4b


4c

4d


4e


4f


Provox Life Home/Go/Energy/Protect HME

5.


6.


7.


Provox Life Night HME

8.


9.





10.

11.


12.

16.

17.

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Contents

EN - ENGLISH LARYTUBE	13
EN - ENGLISH HME	21
DE - DEUTSCH LARYTUBE.....	26
DE - DEUTSCH HME.....	35
NL - NEDERLANDS LARYTUBE ..	41
NL - NEDERLANDS HME	49
FR - FRANÇAIS LARYTUBE	55
FR - FRANÇAIS HME	63
IT - ITALIANO LARYTUBE.....	69
IT - ITALIANO HME.....	77
ES - ESPAÑOL LARYTUBE.....	82
ES - ESPAÑOL HME.....	90
PT - PORTUGUÊS LARYTUBE.....	95
PT - PORTUGUÊS HME	103
SV - SVENSKA LARYTUBE	108
SV - SVENSKA HME	115
DA - DANSK LARYTUBE.....	120
DA - DANSK HME.....	128
NO - NORSK LARYTUBE.....	133
NO - NORSK HME.....	140

FI - SUOMI LARYTUBE.....	145
FI - SUOMI HME.....	153
IS - ÍSLENSKA LARYTUBE	158
IS - ÍSLENSKA HME	165
ET - EESTI LARYTUBE.....	170
ET - EESTI HME.....	178
CS - ČESKY LARYTUBE.....	183
CS - ČESKY HME	191
HU - MAGYAR LARYTUBE	197
HU - MAGYAR HME	205
PL - POLSKI LARYTUBE.....	210
PL - POLSKI HME.....	219
HR - HRVATSKI LARYTUBE	224
HR - HRVATSKI HME	232
EL - ΕΛΛΗΝΙΚΑ LARYTUBE.....	237
EL - ΕΛΛΗΝΙΚΑ HME.....	246
JA - 日本語 LARYTUBE	252
JA - 日本語 HME.....	260
AR - العربية LARYTUBE	265
AR - العربية HME	271

Provox Life System

1. Description of the device

Provox Life devices are designed to work together as a system. Provox Life System consists of Provox Life Adhesives, Provox Life HMEs, Provox Life LaryTubes, Provox Life LaryButtons and Provox Life accessories range of devices.

2. Contraindications

The devices should only be used in accordance with the Instructions for Use.

Patients without the physical, cognitive, or mental ability required to attach, remove or operate the devices themselves, should not use the devices independently and should only use them if they are under sufficient supervision of a clinician or a trained caregiver.

Only valid for HMEs

The devices should not be used by patients with a low tidal volume, as the added dead space may cause CO₂ (Carbon dioxide) retention, see 2.1. Technical data.

Provox Life LaryTube

1. Intended use

Provox Life LaryTube is a single patient use device intended to provide attachment for Provox Life HMEs and accessories after total laryngectomy. For laryngectomized patients with a shrinking tracheostoma

it is also used to maintain the tracheostoma for breathing.

2. CONTRAINDICATIONS

Provox Life LaryTube is not intended to be used by patients that:

- are under any form of mechanical ventilation.
- have damaged tracheal or tracheostoma tissue.

3. Description of the device

Provox Life LaryTube is a tube made of medical grade silicone rubber. It is delivered single packed and non-sterile. Provox Life LaryTube provides attachment for Provox Life HMEs and Provox Life Shower.

There are 4 models of Provox Life LaryTube: (Fig. 1).

Standard model (Fig. 1a) – Made for use with or without a voice prosthesis. Can be attached with a Provox TubeHolder or Provox LaryClips.

Fenestrated model (Fig. 1b) – Made for use with a voice prosthesis. Can be attached with a Provox TubeHolder or Provox LaryClips.

Standard with Ring model (Fig. 1c) – Made for use with or without a voice prosthesis. Can only be attached with a Provox Life Adhesive.

Fenestrated with Ring model (Fig. 1d) – Made for use with a voice prosthesis. Can only be attached with a Provox Life Adhesive.

The Standard and the Ring model can be fenestrated so that air can go through the voice prosthesis. The holes are punched by using Provox Fenestration Punch according to its Instructions for Use. The fenestrations should be in a position that will optimize the airflow through the voice prosthesis.

The different parts of Provox Life LaryTube are (Fig. 1e-f):

e) Fenestration holes

f) Ring for adhesive connection

4. WARNINGS

- DO NOT reuse a used Provox Life LaryTube in another patient. Reuse in another patient may cause cross contamination.
- DO NOT use lubricating gels if the patient has allergies related to these substances.
- DO ensure that the patient has been trained in the use of the device. The patient should demonstrate the ability to understand and consistently follow the Instructions for Use without clinician supervision.

5. PRECAUTIONS

- Always use a Provox Life LaryTube of an appropriate size. Incorrect size may cause tissue damage, bleeding, irritation, obstruction of the voice prosthesis or may cause the stoma to shrink. HME Cassette may exert pressure on the voice prosthesis if Provox Life LaryTube is not sized appropriately.

- DO NOT use force during insertion. It may cause tissue damage, bleeding or irritation and accidental dislodgement of the voice prosthesis (Fig.3) and/or XtraFlange if used. If the product is pushed too far into trachea it may obstruct breathing.
- Fenestrate only with Provox Fenestration Punch. This will ensure that the holes are small and smooth. Fenestration holes that are too large may cause dislodgement of the voice prosthesis or XtraFlange. Fenestration holes with sharp edges may cause formation of granulation tissue.
- DO NOT use a dirty or contaminated Provox Life LaryTube. It may cause infection.
- DO NOT clean the device while inside the stoma.
- DO NOT use the device until it is completely dry (Fig. 4f). Inhalation of disinfectant fumes can cause severe coughing and irritation of the airways.
- DO NOT clean or disinfect by any other method than described as this might cause product damage and patient injury.
- DO NOT use a damaged device since this might cause aspiration of small parts or cause formation of granulation tissue.
- Only use water-soluble lubricants. Oil-based lubricants (e.g. Vaseline) may damage the product and should not be used.
- Stop using the device if the stoma gets irritated or starts to bleed.
- Patients with bleeding disorders or patients on anticoagulants should not use the device if it causes bleeding that re-occurs.

- Stop using the device if granulation tissue forms in the tracheostoma.
- Use only original components. Other accessories may cause product damage and/or malfunction.

6. Instructions for use

6.1 Preparation

Choose correct size

- Choose the right diameter and length (Fig. 2).

Prepare for insertion

- Always inspect the product before each use. Do not use if damaged or dirty.
- Make sure that your hands are clean and dry.
- If needed, lightly lubricate the tube with a small amount of the following recommended water-soluble lubricants: KY Jelly, SurgiLube and Xylocaine (lidocaine).

6.2 Operating instruction

Insertion (Fig. 3):

Provox Life LaryTube Standard and Fenestrated:

Gently insert Provox Life LaryTube into the tracheostoma and attach with Provox LaryClip or Provox TubeHolder (Fig. 1).

Provox Life LaryTube Standard with Ring and Fenestrated with Ring:

Make sure that the adhesive is centered properly so that Provox Life LaryTube does not press against the sides of the stoma. Gently insert Provox Life LaryTube through the adhesive into the tracheostoma (Fig. 1).

NOTE: Insertion and/or removal of Provox Life LaryTube may occasionally cause some slight bleeding, irritation or coughing which usually diminishes over time. Contact your clinician if these symptoms persist.

6.3 Cleaning and disinfection

Clean and disinfect Provox Life LaryTube before each use. Only use drinking water to clean and rinse the device (Fig. 4a-4f).

1. Clean inside of Provox Life LaryTube with Provox TubeBrush or Provox Swab (Fig.4a).
2. The holes of a fenestrated Provox Life LaryTube should be cleaned with a Provox Brush (Fig. 4b).
3. Rinse Provox Life LaryTube, Provox Brush and Provox TubeBrush in water (Fig. 4c).
4. Place Provox Life LaryTube in water-based mild soap for 15 minutes (Fig.4d).
5. Rinse Provox Life LaryTube with water (Fig. 4c).
6. Disinfect Provox Life LaryTube at least once a day (Fig. 4e) with one of the following methods: Ethanol 70% for 10 minutes or Isopropylalcohol 70% for 10 minutes or Hydrogenperoxide 3% for 60 minutes.
7. Let Provox Life LaryTube dry at least for 2 hours or until it looks dry.

7. Storage instructions

Store Provox Life LaryTube in a clean and dry container at room temperature. Protect from direct sunlight.

8. Device lifetime and disposal

Provox Life LaryTube may be used for a maximum period of 6 months. If the device shows any signs of damage it shall be replaced earlier.

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

9. Compatible products

Use only original components. Other products may cause product damage and/or malfunction.

Devices that can be used with Provox Life LaryTube are:

- Provox Life HMEs (Heat and Moisture Exchangers) for patients breathing through a tracheostoma.
- Provox Life Shower is used to prevent water from entering the Provox Life LaryTube during showering.
- Provox TubeHolder, Freevent Neckband or Provox LaryClip is used to keep the Provox Life LaryTube in place.
- Provox TubeBrush, Provox Swab and Provox Brush
- For Provox Life LaryTube with Ring a Provox Life Adhesive is used as attachment and keeping the Provox Life LaryTube in place.

10. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

Provox Life HME

1. Intended use

Provox Life HMEs are single-use heat and moisture exchangers for patients breathing through a tracheostoma.

2. Description of the device

Provox Life HMEs are single-use devices for pulmonary rehabilitation. They come in different levels of humidification, breathing resistance and filtration that makes them suitable for different situations.

The different Provox Life HMEs are:

- Home HME: when taking it easy
- Go HME: when you are out and about
- Energy HME: when physically active
- Protect HME: when you need protection from bacteria, virus, dust and pollen
- Night HME: when sleeping.

2.1 Technical data

	Pressure Drop at 30 l/min*(Pa)	Pressure Drop at 60 l/min*(Pa)	Pressure Drop at 90 l/min*(Pa)	Moisture loss at VT=1000 ml* (mg/l)
Home	60	200	450	19.5
Go	30	100	230	22.5

Energy	15	50	110	23
Protect**	55	180	320	23
Night	65	210	470	18.5

HME internal volume (Dead space)* max 13 ml

* According to ISO 9360

** Bacterial Filtration Efficiency >98%, Virus Filtration Efficiency >98%. Test method adapted from ASTM F2101.

3. Warnings

- Do not exert unintentional physical pressure on the device. Unintentional or accidental closing or compression of the device may obstruct breathing.
- Do not squeeze, or apply excessive force during attachment or when the device is in place. Excessive force can result in the device being pushed through the adhesive coupling, and may obstruct breathing and block airways.
- Do not use non-compatible devices or non-original devices as it may cause personal injury or damage to the device. Do not dismantle or alter the device.
- Do not re-attach a device clogged with mucus, as it may lead to obstructed breathing.
- Do not attach a device that appears slippery by mucus or other fluids, as it will be hard to attach in a correct position.

4. Precautions

- Always test the function of the device prior to use. The lid should immediately return to its open position after releasing it. (Not applicable for Provox Life Night HME).

- Pay attention when the device needs to be removed or attached during the night as the night-time situation could make it more difficult to orient the device in a correct position, and it may also be difficult to detect if contaminated.
- To reduce the risk of water entering the airways, use a Provox Life Shower when showering.
- The device shall not be used beyond 24 hours, is not reusable and shall not be rinsed in water or in any other solutions. This will substantially reduce the function of the device and can increase the risk for infections due to bacterial colonization.
- Do not administer medicated nebulizer treatment over the device since the medication can become deposited in the device.
- Do not use humidifiers or heated humidified oxygen via a mask over the tracheostoma while using the device. The device will become too wet. If oxygen therapy is required, use only non-heated humidified oxygen.

5. Instructions for use

For guidance see;

- Provox Life Home/Go/Energy/Protect HME: pic. 5-7
- Provox Life Night HME: pic. 8-17

5.1 Operating instruction

5.1.1 Provox Life Home/Go/Energy/ Protect HME

To attach

Hold the device and press lightly until the device connects to the attachment device. See fig. 5.

To speak

Press the lid to allow speech. See fig. 6.

To remove

Hold the device and pull sideways until the device disconnects. See fig. 7.

5.1.2 Provox Life Night HME

To attach

Hold the device by the side openings between your thumb and middle finger. With your index finger on the top of the device (where the moon symbol is imprinted), press lightly until the device connects to the attachment device. Carefully rotate the device in the attachment device to verify that it has connected correctly. See fig. 8-15.

To speak

Occlude the side openings simultaneously to allow speech. See fig. 16.

To remove

Hold the device by the side openings between your thumb and middle finger, allowing your index finger to rest on the top of the device (where the moon symbol is imprinted). Pull lightly sideways until the device disconnects. See fig. 8, 9 and 17.

6. Device lifetime and disposal

- HMEs are single use devices that are designed to provide pulmonary rehabilitation both day and night, and are to be replaced in line with daily activities.
- The number and type of HMEs needed during a 24 hour period will vary based on the situations encountered and the frequency of involuntary coughing.

- Typically, one HME is used during the night and about two HMEs during the day. The HME may need to be replaced more often if clogged due to mucus secretions or other contamination. HMEs are not reusable and cannot be rinsed in water or any other solution - this washes out the special salt that is essential for HME effectiveness.
- Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

7. Compatible products

Provox Life HMEs are compatible with Provox Life attachments.

8. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2024 Coloplast A/S. All rights reserved.
Atos Medical AB