PROVOX®

Cleaning and Sterilization
Reinigung und Sterilisierung
Reiniging en sterilisatie
Nettoyage et stérilisation
Pulizia e sterilizzazione
Limpieza y esterilización
Limpeza e esterilização
Rengöring och sterilisering
Rengøring og sterilisering
Rengjøring og sterilisering
Puhdistus ja sterilointi
Puhastamine ja steriliseerimine
Tīrīšana un sterilizācija
Valymas ir sterilizavimas
Čištění a sterilizace
Tisztítás és sterilizálás
Czyszczenie i sterylizacja
Čišćenje i sterilizacija
Καθαρισμός και αποστείρωση
Temizleme ve Sterilizasyon
Очистка и стерилизация
洗淨及び滅菌
세척 및 멸균
清潔和滅菌

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Reprocessing

Reprocessing includes cleaning, disinfection and sterilization. For information regarding assembly and disassembly as well as possible number of repeated reprocessings see the Instructions for Use (IFU) of the specific device.

CAUTION: This manual describes the cleaning, disinfection and sterilization of reusable products manufactured by Atos Medical AB. It is not intended to be used to reprocess single use devices or single patient use devices. Nor is it intended to be used for devices manufactured by another manufacturer than Atos Medical AB. Failure to comply may cause product malfunction and serious patient injury.

Preparation for reprocessing immediately after use

Disassemble the device according to its IFU prior to reprocessing. In order to remove gross debris, the device should be flushed in cold water (10°-20°C/50°-68°F) with e.g. a shower gun in a deep sink immediately after use. Always ensure an optimal flush - too high flow could create infectious aerosols. Water of a temperature higher than 40°C/104°F should not be used since it could disable the removal of proteins. The device should be cleaned at the reprocessing facility within 24 hours after use.

Precaution

Only use recommended detergents/disinfectants. Some detergents/disinfectants may weaken, damage or destroy the product while others may be absorbed by the product causing the patient to experience adverse tissue reactions. Reprocessing by other methods and with other detergents/disinfectants than those recommended below is under the sole discretion of the user. Always ensure that the disinfectant chosen does not disable further removal of proteins.

Cleaning and disinfection

General

The purpose of the cleaning is to remove debris from the device. Cleaning and disinfection can be done both manually and in an automatic washer-disinfector. The device should be positioned in the washer-disinfector as to achieve the most efficient flushing. This also applies for manual cleaning.

Cleaning validation has shown that the drinking water flow during rinsing should be at least 2 liters/minute.

The device has been validated with and tested to withstand the following alkaline detergent containing sequestering agent:

- Suma Med Super LpH® (Johnson Diversey) and the following enzymatic detergent:
- 3E-Zyme® (Medisafe)

The purpose of the disinfection is to reduce the number of microorganisms on the device.

The device has been validated with and tested to withstand the following chemical disinfectant:

- Gigasept® Instru AF (Schülke&Mayr)

Validated cleaning and disinfection procedures are described below:

Cleaning and disinfection procedure

The device can be cleaned and disinfected by using a washer-disinfector (A) or by manual cleaning and disinfection (B). Disinfection can be performed either chemically or thermally.

A. Washer disinfector

1. Flush the device with cold (10°-20°C)/ (50°-68°F) drinking water for 1 minute.
2. Flush the device with 25°C/77°F drinking water for 1 minute.
3. Clean the device for 10 minutes at the temperature and concentration recommended by the cleaning agent manufacturer.
4. Rinse the device thoroughly with 35°C/95°F drinking water for 1 minute. Repeat the rinsing with drinking water twice (3 rinsing's in all).

5. **Disinfection**

Disinfection can be performed either chemically or thermally.

5.1 Chemical disinfection

1. Disinfect the device chemically at the concentration, temperature and time recommended by the manufacturer of the disinfectant.

2. Rinse the device thoroughly with 35°C/95°F drinking water for 1 minute. Repeat the rinsing with drinking water once. Perform a final rinse for 1 minute with high-purity water, e.g. distilled water.

5.2 Thermal disinfection

Rinse the device with water of a temperature of at least 90°-95°C/194°-203°F for 1 minute. Ensure that the rinse is performed with high purity water e.g distilled water.

6. Let the device dry at 75°C/167°F for 30 minutes. If necessary, continue drying at room temperature in an area with no risk of contamination until the device is completely dry.
7. Inspect the device with regards to residues left after cleaning and disinfection visually or e.g. in a microscope. Lumens could be tested with a fast protein test like the ninhydrin and pyromol test. If the device is not visually clean, repeat cleaning and disinfection. Ensure mechanical integrity of the device before sterilization.

B. Manual cleaning and disinfection

1. Pre-soak the device in cold (10°-20°C)/(50°-68°F) drinking water for 10-20 minutes. Remove debris by flushing the device with 25°C/77°F drinking water for 1 minute.
2. Immerse the device in a cleaning solution at the temperature and concentration recommended by the cleaning agent manufacturer for 10-20 minutes. Clean thereafter mechanically for 5-10 minutes

by using a soft bristled brush and/or a sponge. All surfaces should be cleaned mechanically.

3. Rinse the device with 35°C/95°F drinking water at a flow rate of at least 2 liters/minute for 3 minutes.

4. **Disinfection**

Disinfection can be performed either chemically or thermally.

- 4.1 **Chemical disinfection**

1. Disinfect the device by immersing it in a chemical disinfectant solution at the concentration, temperature and time recommended by the manufacturer.

2. Rinse the device with 35°C/95°F drinking water at a flow rate of at least 2 liters/minute for 2 minutes. Repeat the rinse for 1 minute with high-purity water, e.g. distilled water.

- 4.2 **Thermal disinfection**

Immerse the device in water that is 90°-95°C/194°-203°F for 1 minute. Ensure that the disinfection is performed with high-purity water e.g. distilled water.

5. Let the device dry at room temperature in an area with no risk of contamination until the device is completely dry.
6. Inspect the device with regards to residues left after cleaning and disinfection visually or e.g. in a microscope. Lumens could be tested with a fast protein test like the ninhydrin and pyromol test. If the device is not visually clean, repeat cleaning and disinfection. Ensure mechanical integrity of the device before sterilization.

Sterilization

The device has been tested to withstand steam sterilization. Make sure that the device has been cleaned, disinfected and is completely dry prior to sterilization. Ensure that the device is assembled/disassembled according to the IFU of the device. Put the device in sterilization pouches permeable for steam and air.

Dynamic air-removal steam sterilization procedures have been validated with the overkill method described in EN ISO 17665-1. Accordingly, the recommended sterilization times below are twice the sterilization times used in the validation.

Sterilization at 121°C/249.8°F for 25 minutes after performing three 4 minutes vacuum pulses with depths of 0.06 bar.

Sterilization at 134°C/273.2°F for 3 minutes after performing three 4 minutes vacuum pulses with depths of 0.06 bar.

In order to achieve dry devices, drying after sterilization was performed for 8 minutes at a temperature of at least 80°C/176°F during a pressure reduction down to 0.04 bar.

The steam used was of grade according to SS-EN 285.

Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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© Coloplast A/S. All rights reserved.
Atos and the Atos Medical logo are trademarks of Coloplast A/S.